



MATU INTERNATIONAL PVT. LTD.

B-11 Jhilmil Industrial Area Delhi-110095,INDIA

Customer care no. 9821133428

Email-sales@matuhealthcare.com



BLOOD TRANSFUSION SET

Features :

- o For intra vascular, gravity controlled transfusion of blood or blood elements.
- o Soft, flexible and transparent drip chamber.
- o Soft and kink resistant DEHP free P.V.C. tube.
- o For approximately 20drops/ml.
- o With 200 micron blood filter.

Option Available:

- o With or Without Needle.
- o Luer Lock or Luer Slip.
- o With airvent spike/without airvent spike..
- o Latex Bulb for intermi
- o Single/Double Drip chamber.
- o Packed in pouch or blister

TECHNICAL SPECIFICATION :

SIZE	CAPACITY	PACKING INNER/OUTER
VENTED	20 Drop/ml	25/500
NON-VENTED	20 Drop/ml	25/500