











MATU INTERNATIONAL PVT. LTD.

B-11 Jhilmil Industrial Area Delhi-110095,INDIA
Mobile no. 9821133428 | Email-sales@matuhealthcare.com

Contact : Nishant BD manager Telephone/Whatsapp: +91-8130573898 sales@matuhealthcare.com	OEM Product List	Trade terms: FOB INDIA ICD TKD New Delhi	CE ISO13485
Name	Infusion Sets IV	Picture	Product Code
Infusion sets for single use	1. LUER SLIP: Export basic type,48mm chamber with air vent with fluid filter,1.5M standard PVC tube(2.9*3.75mm ID/OD),Standard PP/PE flow regulator, rubber bulb, LuerSlip connector with Hypodermic Needle		A1
Infusion sets for single use	2. LUER LOCK: Export basic type,48mm chamber with air vent with fluid filter,1.5M standard PVC tube(2.9*3.75mm ID/OD), Standard PP/PE flow regulator, rubber bulb, LuerLock connector with Hypodermic Needle		A2
Infusion sets for single use	3. LUER SLIP: Export basic type,48mm chamber with air vent with fluid filter,1.5M standard PVC tube(2.9*3.75mm ID/OD),Standard PP/PE flow regulator, rubber bulb, LuerSlip connector with Protetive Cap		A3
Infusion sets for single use	4. LUER LOCK: Export basic type,48mm chamber with air vent with fluid filter,1.5M standard PVC tube(2.9*3.75mm ID/OD),Standard PP/PE flow regulator, rubber bulb, LuerLock connector with Protetive Cap		A4

Contact : Nishant BD manager Telephone/Whatsapp: +91-8130573898 sales@matuhealthcare.com	OEM Product List	Trade terms: FOB INDIA ICD TKD New Delhi	CE ISO13485
Name	Infusion IV Sets (Y-Site)	Picture	Product Code
Infusion sets for single use	1. LUER SLIP: Export Y - type,48mm chamber with air vent with fluid filter,1.5M standar d PV C tube(2.9*3.75mm ID/OD) , Standar d PP/PE flow regulator, Y -shaped liquid medicine adding device with Connecting catheter, LuerSlip connector with Hypodermic Needle		Y1
Infusion sets for single use	1.LUER LCOK: Export Y - type,48mm chamber with air vent with fluid filter,1.5M standar d PV C tube(2.9*3.75mm ID/OD) , Standar d PP/PE flow regulator, Y -shaped liquid medicine adding device with Connecting catheter, LuerLock connector with Hypodermic Needle		Y2
Infusion sets for single use	1. LUER SLIP: Export Y - type,48mm chamber with air vent with fluid filter,1.5M standar d PV C tube(2.9*3.75mm ID/OD) , Standar d PP/PE flow regulator, Y -shaped liquid medicine adding device with Connecting catheter, LuerSlip connector with Protetive Cap		Y3
Infusion sets for single use	1.LUER LCOK: Export Y - type,48mm chamber with air vent with fluid filter,1.5M standar d PV C tube(2.9*3.75mm ID/OD) , Standar d PP/PE flow regulator, Y -shaped liquid medicine adding device with Connecting catheter, LuerLock connector with Protetive Cap		Y4

Contact : Nishant BD manager Telephone/Whatsapp: +91-8130573898 sales@matuhealthcare.com	OEM Product List	Trade terms: FOB INDIA ICD TKD New Delhi	CE ISO13485
Name	Infusion IV Sets (Y-Site)	Picture	Product Code
Infusion sets for single use	Export Y - type, 48mm chamber with air vent with fluid filter, 1.5M standard PV C tube (2.9*3.75mm ID/OD), Standard PP/PE flow regulator, Y-shaped liquid medicine adding device with Connecting catheter, LuerSlip connector with butterfly Scalp Vein Set		Y5
Infusion sets for single use	Export Y - type, 48mm chamber with air vent with fluid filter, 1.5M standard PV C tube (2.9*3.75mm ID/OD), Horn shape/Standard PP/PE flow regulator, Y-shaped liquid medicine adding device with Connecting catheter, LuerLock connector with butterfly Scalp Vein Set		Y6
Name	Infusion IV Sets (Classic-Style)	Picture	Product Code
Infusion sets for single use	Classic style + New air filter device. ABS Plastic apike needle + New Air filtration Set + chamber + 1.5M standard PVC tube (2.9*3.75mm ID/OD) + Standard PP/PE flow regulator + LuerLock connector with fluid filter + Scalp Vein Set		C1



MATU INTERNATIONAL PVT. LTD.

B-11 Jhilmil Industrial Area Delhi-110095,INDIA

Mobile no. 9821133428

Email-sales@matuhealthcare.com

Contact : Nishant BD manager Telephone/Whatsapp: +91-8130573898 sales@matuhealthcare.com	OEM Product List	Trade terms: FOB INDIA ICD TKD New Delhi	CE ISO13485
Name	Infusion IV Sets (Classic-Style)	Picture	Product Code
Infusion sets for single use	Classic style+New air filter device. ABS Plastic apike needle + New Air filtration Set+ chamber + 1.5M standard PVC tube (2.9*3.75mm ID/OD)+ Standard PP/PE flow regulator+ LuerLock connector with fluid filter + butterfly Scalp Vein Set		C2
	Otherwise: 60 drops/ml infusion set for children. Various customized infusion sets, please send pictures for inquiry.		

Explanation:

- I Support customization.
- 1 According to actual demand,each part of the infusion set product can provide customized combination scheme and provide corresponding quotation.
- 2 Packaging also supports OEM customization.
- 3 Product validity: 5 years.
- II About supply capacity: Each specification can supply 12 million Pcs per month.
- III MOQ:450000 PCS.
- IV Carton Size:62*38*33cm,25Pcs/bag,500Pcs/CTN
(A 40HQ Container: Hold 450000 PCS)

V Basic commercial terms:

- 1 Under FOB terms: The prices is FOB INDIA ICD TKD NEW DELHI
- 2 Place of delivery : ICD TKD NEW DELHI,INDIA
- 3 The personnel designated by the buyer shall inspect in the factory warehouse, and deliver the goods after confirmation. Before each batch of goods leaves the factory warehouse, the final payment should be settled.