



Manufactured By: **MATU INTERNATIONAL PVT. LTD.**

B-11 Jhilmil Industrial Area Delhi-110095

Customer care no. 9821133428

Email-sales@matuhealthcare.com

**MADE IN INDIA**

**People-oriented, technological innovation, harmonious development**



Color Codes and Specifications

O, D (mm)	Gauge	Color Code			13mm	16mm	19mm	25mm	32mm	38mm
					1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"
0.36	28G	Blue-green								
0.40	27G	grey								
0.45	26G	brown								
0.50	25G	orange								
0.55	24G	purple								
0.60	23G	blue								
0.70	22G	black								
0.80	21G	green								
0.90	20G	yellow								
1.00	19G	cream								
1.20	18G	pink								
1.60	16G	white								



## Syringes



Sterile Syringe for Single Use



2-part Syringe



3-part, Luer Lock



Needle